

# Peak Neptune Self-Propelled Wheelchair

Product Code: 9346376060516 (Red)













## Features:

- ▶ Red Aluminum Frame
- ▶ Drop-Back Handle
- ▶ Flip-Up Armrest
- ▶ Detachable Footrest
- ▶ Durable PU Rear Wheel
- ▶ Easily Folds



## ARTG Listing no.

- ▶ 313946

Length (cm)	 100	Backrest Height (cm)	 42	Front Wheel Diameter (cm)	 20/8"
Width (cm)	 61.5	Seat Depth (cm)	 40	Rear Wheel Diameter (cm)	 60/23"
Height (cm)	 92	Seat Width (cm)	 46	Max Load (kg)	 100
Folded Width (cm)	 28	Seat Height from the Ground (cm)	 47	Net Weight (kg)	 13.1