

Peak Care Electric Home Care Bed Pull Apart

Product Specifications 2016

Product Codes: BH-989 & BH-989MABS



Specifications

- Australian Standards Tested to AS/NZS 3200.1.0:1998 & AS/NZS 3200.2.38:2007 to 200kg
- Weight capacity 200kg
- Overall size: L2300 x W1000 mm
- Knee Brake Adjustment 0° – 22°
- Adjustable Bed Height 390 – 790 mm
- Trendelenburg +/- 10°
- Steel frame made with special antibacterial powder coating
- Steel Plate Sleep Surface
- Low Voltage DC T-Motion Motor System
- 4 x High Quality Lockable Castors
- Motor Covers Standard
- Code BH-989MABS comes with ABS Side Rails

ARTG Listed
No. 196516

Tested by a NATA
accredited laboratory
Accreditation No:
2953



Quality Assurance